



# ACCESSORIES

## I. V. POLES

For infusion pumps and fluids bags



### M-IS16

- For 4 infusion fluids bags/bottles
- Up to 6 infusion pumps
- Stainless steel construction
- Rail-mounted
- Height adjustable
- Various length of turning arms up to 500mm
- Maximum load - 20kg.

### Sliding I.V. pole for Pendant Bridges M-GM

- For 4 infusion fluid bags/bottles
- Stainless steel construction
- On special linear rail system with movement brakes and additional railings for extra accessories
- Height adjustable
- Various length turning arm up to 500 mm
- Maximum load - 25kg.



I.V. Holder for Bed Head Unit M-PB



Infusion pump pole  
for rail mounting



## RAILINGS FOR MOUNTING DEVICES

### M-RL05

- Wall mounted rails for suspended equipment
- Standard rail dimensions 25 x 10 mm, 30 x 10 mm
- Maximum load 40kg/meter
- Stainless steel construction



- For mounting to Bed Head Unit profile



- For various suspended equipment

## BEDSIDE SCREENS



M-BS05T

- Two types of holders:
  - telescopic M-BS05T
  - ceiling profile M-BS05C
- Various length and height holders
- Large selection of curtain materials
- Easy cleaning
- Flame resistant materials
- Fixation of telescopic bedside screens:
  - on the wall
  - on the rails
- Special profile for suspended ceiling



M-BS05C

## ACCESSORIES FOR MAIN EQUIPMENT



### Sliding shelf module for Pendant Bridges M-GM

- Mobile system can be moved along the pendant with rotation about vertical axis, with a separate height-adjustable shelf
- On special linear rail system with movement brakes
- Made from stainless steel and aluminum, which is painted with powder paint in metallic colors
- Dimensions may be standard or by customer's request;
- Resistant to all the cleaning detergents
- Maximum load - 40kg/per self

### Suspended shelf with drawers for Wall Pendants M-GB

- Aluminum painted with powder paint with the top shelf made from stainless steel
- Easy to clean, without sharp edges
- Operative mounting to the profile of wall pendant
- Resistant to all the cleaning detergents
- Maximum load - 40kg.



### Shelf with drawer for Ceiling Pendants M-GL



- Two rails 10x25 for suspended equipment
- Handle to control pendant's pneumobrakes
- Made from aluminum, which is painted with powder paint in metallic colors
- Dimensions may be standard or by customer's request
- Resistant to all the cleaning detergents
- Maximum load - 40kg.

### Shelf for Ceiling Pendants M-GL

- Rail 10x25 for suspended equipment
- Possibility to adjust the height of mounting
- Made from stainless steel and aluminum, which is painted with powder paint in metallic colors
- Dimensions may be standard or by customer's request
- Resistant to all the cleaning detergents
- Maximum load - 40kg.



### Monitor holder



- For various sizes monitors fixation to the ceiling pendants, wall pendants, bed head units.
- Maximum load 17 kg
- Arm length up to 500 mm



### Shelf for the monitor for Wall Pendants M-GB



- Comfortable mounting of the shelf through the whole length of the pendant
- Height adjustment
- Mounting on the rail 10x25 mm
- Made from aluminum, which is painted with powder paint in metallic colors
- Dimensions may be standard or by customer's request
- Resistant to all the cleaning detergents
- Maximum load - 15kg.

### Shelf for the monitor for Bed Head Unit M-PB

- Comfortable mounting of the shelf through the whole length of the pendant
- Made from stainless steel and aluminum, which is painted with powder paint in metallic colors
- Dimensions may be standard or by customer's request
- Resistant to all the cleaning detergents
- Maximum load - 15kg.



Email: [mtlbi@medtech.lt](mailto:mtlbi@medtech.lt)