



MULTIPOINT™ TERMINAL UNITS

PENDANT BRIDGES M-GM

FEATURES

- Supply of medical gases, compressed air and vacuum
- Supply of electricity, low-voltage and audio-video signal
- Easy fitting of control and life support devices on shelves, rails and I.V. poles
- Possibility of installation in any desired location in the facility, like in front of windows, glass partitions
- Easy to access patient from all sides
- Prompt fitting of equipment on the sliding shelf module and I.V. pole.
- Special four-sections anodized aluminum profile 250 x 116 mm, length depends on a project
- The front panels and side covers are from anodized aluminum or coated with powder paint in RAL colors
- Two railings 10 x 25 mm for accessories
- Indirect lighting, patient's zone lighting and night light
- Medical gas outlets, electric sockets, information signal sockets, medical gas pressure control indicators





SET OPTIONS:

Pendant Bridge can be equipped with:

- Sliding shelf module can be moved along the pendant with rotation about vertical axis with a separate height-adjustable shelf
- A shelf for a monitor
- I. V. pole for infusion pumps and drips, with standard rail 10x25 mm for suspended equipment
- Fixation for fluid collection jars
- Sources of direct light for patient's examination
- Bedside telescopic screens
- Flowmeter with oxygen humidifier mounted on a rail or with direct connection
- Systems of suction from centralized vacuum, compressed air with fluid collection jars
- Catheter holder and other suspended devices
- Remote control for nurse call and lighting system control

FIELD OF APPLICATION



- For small operating wards
- For narcosis wards
- For intensive therapy wards and resuscitation rooms
- For neonatal wards
- For emergency rooms
- For treatment rooms



Pendant Bridges can be joined end-to-end to form a continuous module for general resuscitation rooms, intensive therapy wards and others.

CONSTRUCTION FEATURES

- Solid fixation of profile onto the ceiling with two support-legs
- Four-section special anodized aluminum profile with two special rails for accessories
- Accesories can be moved along the pendant:
 - sliding shelf module
 - sliding I.V. Pole
 - shelf for monitor
- Front covers are made from aluminum profile and fixed into main profile without any fasteners
- Side covers made with pins for internal fixation



- Profile length from 1500 mm up to 2800 mm according to the project
- Separate medical gas and electric sections for connection with inlets through separate sections

Colors of front panels and side covers:



Color of profile: Color of shelves:



2 YEARS GUARANTEE for all pendants bridges.

Pendants bridges are designed and manufactured in accordance with directive 93/42/EEC requirements.

Email: mtlbi@medtech.lt