



MOBILE PROCEDURAL

FEATURES



800 600 1050 mm

- 100% polished stainless steel
- Removable shelves with edges
- Casters with brakes
- Two fixed positions drawer
- Dismountable, easy transportable

TECHNICAL DATA

- Dimensions of shelves
750 570 15 mm
650 470 15 mm
- Load capacity of one shelf up to 20 kg



700 500 1050 mm

MOBILE INSTRUMENT TABLE M-IS12

FEATURES



- 100% polished stainless steel
- Height-adjustable shelf
- Easy removable tray with edges
- Increased weight of chassis to ensures stability of table
- Casters with brakes
- Fixed position of table relative to the chassis
- Dismountable, easy transportable

TECHNICAL DATA

- Two sizes of tray table:
500 x 400 mm
700 x 400 mm
- Table height adjustment from the floor 900 ? 1300 mm
- Load capacity of tray up to 20 kg

WALL MOUNTED BEDSIDE CABINETS FOR PATIENTS



TECHNICAL DATA

- Two different models:
- M-MPC03 - standard cabinet for storage of personal belongings
Dimensions - 605 460 372 mm
Dimensions of the drawers:
369 322 180 mm.
396 322 252 mm
- M-MPC05 - with adjustable over-bed table
Dimensions - 605 460 372 mm
Dimensions of the drawers:
369 322 110 mm.
396 322 252 mm
Over-bed table adjustment range:
Lifting height of over-bed table from 750mm to 945mm (when cabinet is mounted 150mm over floor)
Maximum radius of the over-bed table from wall - 1075mm

FEATURES

- Suspended above the floor for easy cleaning
- For various type of wards
- Comfortable feeding of serious ill patients using over-bed table
- Two different height drawers
- Shelf for shoe storage



1 YEAR GUARANTEE for all types of tables

Email: mtlbi@medtech.lt