





CENTRAL GAS SUPPLY

OXYGEN SUPPLY PLANT

- Oxygen supply station with liquid oxygen from cryogenic vessels as a primary source and two sources from oxygen collectors
- Cost effective in hospitals with 100 beds and more



- As a primary source of oxygen supply in hospitals can be used oxygen concentrator instead.
- Oxygen supply from cylinder source may be used for supplying oxygen to hospitals up to 50 beds or



HIGH PRESSURE COLLECTORS

- With fixation of cylinders to the wall from 2x2 to 2 x10 cylinders
- M-R2N - for nitrous oxide
- M-R2C - for carbon dioxide
- M-R2O - for oxygen




- With fixation of cylinders to the special frame from 2x8 and 4x8 cylinders
- M-R3N - for nitrous oxide
- M-R3C - for carbon oxide
- M-R3O - for oxygen

GAS PRESSURE REDUCTION AND CONTROL UNITS

RPA 1500 K

- For connection of three supply sources
- One inlet 7 - 20 bar
- Two inlets 10 - 200 bar
- Outlet pressure 3 - 10 bar
- Two parallel branches for pressure reduction



- Productivity up to 1500 l/min
- Digital pressure control
- Light and sound alarm system
- Heated high pressure reducers
- Safety valves

RPA 1500 B



- Two high pressure inlets 10 - 200 bar
- Two stage pressure adjustment
- Outlet pressure 3 - 10 bar
- Productivity up to 1500 l/min
- Light and sound alarm system
- Heated high pressure reducers
- Safety valves

RP-500




- For one manifolds
- Two stage pressure adjustment
- Heated high pressure reducer
- Safety valve
- Pressure control of supply system
- Productivity up to 500 l/min

RPA-500S

- For two collectors with automatic changover
- Two stage pressure adjustment
- Heated high pressure reducers
- Safety valves
- Pressure control, transmission of alarm signal
- Indication of supply source
- Productivity to 500 l/min



TM



- Semi-flexible stainless steel tubes with various types brass joints for connection of collectors and cylinders
- Maximum pressure up to 30 Mpa
- Connection joints for cylinders according to client's national standard

- Semi-flexible thick-walled copper tubes for connection of reducing device with high pressure collectors
- Maximum pressure to 30 Mpa



1 YEAR GUARANTEE for all medical gas cylinder stations

TM

Email: mtlbi@medtech.lt